

# Caring for Life's Greatest Treasures



Introducing our Pediatric Blanket and Fluid Warmers  
Make sure to request the label!



\*Blanket Warmers supplied with dog/cat label

\*\*Fluid Warmers supplied with goldfish label

Compact/Mobile Blanket Warmers



EC/DC250

- Capacity: 4-6 blankets
- 22" H; 18.5" W; 23.5" D
- Available via UPS shipment
- Countertop model
- Available with casters
- Optional cart

EC/DC350

- Capacity: 10-12 blankets
- 28" H; 18.5" W; 23.5" D
- Floor or countertop model
- Available with casters
- Optional cart

EC/DC400

- Capacity: 14-18 blankets
- 22" H; 24" W; 25.75" D
- Floor or countertop model
- Available with casters
- Optional cart

EC/DC750

- Capacity: 14-18 blankets
- 35" H; 24" W; 25.75" D
- Floor or countertop model
- Available with casters

\*Noted depth for EC Series units including the handle

Fluid Warmers - Cleared for IV and irrigation fluids



EC/DC250L

- Capacity: 16 (1) liter bags or 12 (1) liter bottles
- 22" H; 18.5" W; 23.5" D
- Available via UPS shipment
- 510k cleared
- Countertop model
- Available with casters
- Optional cart

EC350L

- Capacity: 24 (1) liter bags or bottles
- 28" H; 18.5" W; 23.5" D
- 510k cleared
- Countertop model
- Available with casters
- Optional cart

EC/DC400L

- Capacity: 28 (1) liter bags or 20 (1) liter bottles
- 22" H; 24" W; 25.75" D
- 510k cleared
- Countertop model
- Available with casters
- Optional cart



MS-400

- Fits 250, 350 and 400 size models
- 36.5"H ; 26.0"W; 26.0"D