

TITAN SERIES

Affordable, precise temperature control

HEAVY DUTY WARMERS

EC250L

- Total capacity: 12 1-liter bottles or 16 1-liter bags
- Exterior size: 22" H; 18.5" W; 23.5" D
- Countertop model
- Stackable



EC400L

- Total capacity: 20 1-liter bottles or 28 1-liter bags
- Exterior size: 22" H; 24" W; 25.75" D
- Countertop model
- Stackable



EC350L

- Total capacity: 24 1-liter bottles or 24 1-liter bags
- Exterior size: 28" H; 18.5" W; 23.5" D
- Floor or countertop model
- Stackable



EC770L


- Total capacity: 54 1-liter bottles or 72 1-liter bags
- Exterior size: 41" H; 25.25" W; 35.5" D
- Stackable



ivNext®

iv Inventory Management

INTELLIGENT FLUID WARMER



Ready to Use

Warming

Expired

- Total capacity: 10 1-liter bags
- Exterior size: 22" H; 18.5" W; 23.5" D
- Floor or countertop model
- Stackable



Blue = Warming



Green = Ready to Use



Yellow = Expired

Organized inventory allows for faster retrieval of fluids

- First-in, First-out
- Know when an iv bag has reached set point temperature
- Visual indication of bag readiness or expiration
- Quickly retrieve temperature documentation with standard data logging feature

Simple for the clinician
Improve Patient Safety
Increase regulatory compliance

Increase regulatory compliance
Save money
Improve processes
Increase reimbursements

ivNow®

FLUID WARMERS



The unit shown above provides 75 liters of warm fluid in a 12 hour period.

Warmed fluid where you need it, when you need it

- Know the time and temperature of every IV bag
- Warms 1L iv bags in 30 minutes or less
- Warms 25-78 liters in 12 hours
- Point-of-care use saves space, time & footsteps
- Requires no disposables



COMFORT SERIES

WARMS SAFELY, QUIETLY & EVENLY

DC250L

- Total capacity: 12 1-liter bottles or 16 1-liter bags
- Exterior size: 22" H; 18.5" W; 22" D
- Countertop model
- Stackable



DC400L

- Total capacity: 20 1-liter bottles or 28 1-liter bags
- Exterior size: 22" H; 24" W; 24.3" D
- Floor or countertop model
- Stackable



MOBILE CARTS

MS-400

- Mobile cart for small blanket and fluid warmers
- Fits the following models: DC250, DC250L, DC350, DC400, DC400L, EC250, EC250L, EC350, EC350L, EC400, EC400L, ivNext
- Exterior size: 36.5" H; 26" W; 26" D



Shown with DC250L Fluid Warmer

MS-550

- Mobile cart for Multi-Purpose Warmers
- Fits the following models: EC550BL-18" Depth; EC550BL-24" Depth
- Exterior size: 34" H; 32.5" W; 29.5" D



Shown with EC550BL 18" Depth Combination Warmer

FLUID WARMERS

YOU CAN AFFORD THE BEST



enthermics®
An ISO 13485:2016 certified company

www.enthermics.com

RETHINK PATIENT WARMING



More Than Just Comfort

Patient warming is about more than providing comfort. Maintain normothermia throughout a procedure and hospital or clinical stay to:

- Decrease blood loss
- Reduce incidence of surgical site infection
- Shorten stays and reduce readmissions
- Improve patient satisfaction and comfort

Enthermics' flexible warming technologies help you maintain patient normothermia from admission to discharge. Available in a range of sizes and types, we have warming solutions to fit most any budget.



ivNow-1



EC250L over EC350

MULTI-PURPOSE WARMER

EC550BL-18

- Total capacity:
 - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 9-12 blankets
- Exterior size: 25.12" H; 30.0" W; 21.12" D
- 18" chamber depth
- Also available as a built-in cabinet

Unit shown with fluid bottles, bags and blankets for illustration purposes only



EC550BL-24

- Total capacity:
 - TriMode chamber: 30 1-liter bottles + 36 0.5-liter bottles, 48 1-liter bags or 10-14 blankets
- Exterior size: 25.12" H; 30.0" W; 28.5" D
- 24" chamber depth
- Also available as a built-in cabinet

Unit shown with optional basket - set up for fluids



COMBINATION BLANKET & FLUID WARMERS (DUAL CAVITIES)

EC1350BL

- Total capacity:
 - Blanket chamber: 24-30 blankets; 9.39 ft³ interior size
 - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 9-12 blankets
- Exterior size: 75" H; 30.0" W; 21.25" D
- 18" chamber depth
- Also available as a built-in or a pass-through cabinet
- OSHPD-OPM-0562-19



Shown as optional pass-through cabinet

EC1850BL

- Total capacity:
 - Blanket chamber: 30-40 blankets; 13ft³ interior size
 - TriMode chamber: 30 1-liter bottles, 24 1-liter bags or 10-14 blankets
- Exterior size: 75" H; 30.0" W; 28.5" D



EC1540BL

- Total capacity:
 - Blanket chamber: 20-25 blankets; 7.7 ft³ interior size
 - Fluid chamber: 54 1-liter bottles or 72 1-liter bags
- Exterior size: 72.75" H; 24" W; 32.5" D



STORAGE CABINET

1560-GS

- Interior size: 15.6 ft³
- Exterior size: 75" H; 30" W; 21.25" D
- Pairs with EC1350BL Combination Warmer




MIX & MATCH

Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures

Match the right temperature to the right product in the right place at the right time.

Only Enthermics offers such a full line of fluid and blanket warmers to stack to create the warming configuration to fit your needs.



More combinations available

enthermics



5-Year Parts, 1-Year Labor Warranty on cabinets
* US/Canada only

ENTHERMICS MEDICAL SYSTEMS

An ISO 13485:2016 certified company

Tel 262-251-8356 | 800-TO-B-WARM

generalinfo@enthermics.com

www.enthermics.com