

Carts



The MS-400 and MS-550 carts are mobile carts designed to hold the Enthermics table-top warming cabinets.

Convenient: Mobile

Service: Our worldwide support team is dedicated to providing quality equipment. We walk with you every step of the way to maximize your return on investment.

MS-400 fits the following models: DC250, DC250L, DC350, DC400, DC400L, EC250, EC250L, EC350, EC350L, EC400, EC400L, ivNext®.

MS-550 fits the following models: EC550BL-18" depth, EC550BL-24" depth

Made in the USA with a commitment to quality

An ISO 13485:2016 certified company

Standard Features

- Stainless Steel
- 3-1.2" (89mm) casters, 2 rigid, 2 swivel with brake
- 1 shelf



MS-400 Cart with EC250L







MS-550 Cart with EC550BL-18" Depth

EC550BL shown with iv bags, irrigation bottles and blankets for illustration purposes only. Warming mode should be selected based on content!.

MS-400

MS-550



Model	Exterior (H x W x D)
MS-400	36.5" (927mm) x 26" (660mm) x 26" (660mm)
MS-550	34" (866mm) x 32.5" (825mm) x 29.5" (749mm)

Net Weight
50 lbs (23 kg)
56 lbs (25 kg)

	Ship Dimensions (L x W x H)*	Ship Weight*
MS-400	35" x 35" x 44"	120 lbs
MS-550	35" x 35" x 44"	125 lbs

*Domestic ground shipping information. Contact factory for export weight and dimensions.

CONTACT US

W164 N9221 Water Street | Menomonee Falls, Wisconsin 53051 | U.S.A.
 Phone: 262.251.8356 | 800.862.9276 | Fax: 262.251.7067 | enthermics.com