

Streamline Your Processes

The Benefits of Point-of-Care Fluid and Blanket Warmers



Infection Prevention

- Having warmed fluids and blankets at point-of-care reduces exposure to harmful viruses and other contaminants.

Save time and money

- Fluids and blankets available at point-of-care.
- No special set up costs.

Improve efficiencies/Save footsteps

- No need to leave the room and go to other departments to get warm fluids or blankets.
- Not having to leave the O.R. can ensure that the sterile environment is not compromised.
- Patients are assured warm fluids and blankets as there is no heat loss due to transport.

Improve compliance/JAHCO, AAAA etc. by utilizing the ivNow or ivNext. *The joint commission requires that fluid warmers are monitored daily for temperature and inventory management.*

- The ivNow® precisely displays the actual temperature while monitoring the shelf life of each bag while the ivNext® indicates bag readiness via color-coded lighting.

Improve Clinical Outcomes/SSI's (surgical site infections)

- Providing your patients with warmed fluids reduces the risk for peri-operative hypothermia which reduces their risk for SSI's by up to 66%.




Fluid Warmers - Cleared for IV and/or irrigation fluids



EC250L

- Capacity: 16 1-liter bags or 12 1-liter bottles
- 22" H; 18.5" W; 23.5" D
- 510K cleared
- Stainless Steel
- 1 basket standard
- Countertop model
- Stackable



ivNext®

- Capacity: 10 1-liter bags
- 22" H; 18.5" W; 23.5" D
- 510K cleared
- Stainless Steel
- 10 shelves
- Countertop model
- Stackable



DC250L

- Capacity: 16 1-liter bags or 12 1-liter bottles
- 22" H; 18.5" W; 22" D
- 510K cleared
- Painted Steel
- 2 baskets standard
- Countertop model
- Stackable

ivNow®

- Capacity: Will provide up to 75 liters of warmed fluid every 12 hours
- No training or disposables required
- 510K cleared
- Easy-to-Use: Place the IV bag in the pod to activate the unit. Fluid is automatically warmed up to 104°F (40°C) in 30 minutes or less.
- Configurations: Countertop, mobile equipment pole stand, boom or wall.



ivNow-1



ivNow-2



ivNow-3

Compact/Mobile Blanket Warmers



EC250

- Capacity: 4-6 blankets
- 22" H; 18.5" W; 23.5" D
- Stainless Steel
- Countertop model
- Stackable



DC250

- Capacity: 4-6 blankets
- 22" H; 18.5" W; 23.5" D
- Painted Steel
- Countertop model
- Stackable