

Titan Series

Storage Cabinet

Store your supplies for busy emergency rooms, operating rooms and labor and delivery units in the convenient Titan Series storage cabinet

The durable, heavy-duty stainless steel Titan Series storage cabinet allows you to store your blankets and fluids next to your warming cabinet to easily replenish the warming cabinet as needed. Partners directly with the EC1350BL Combination Warmer

Durable: The affordable, rugged stainless steel cabinets are made to last and work just as hard as you do.

Convenient: Easily replenish your warming cabinet with supplies from the storage cabinet.

Glass Window: The glass window allows staff to take an immediate inventory of the amount of available supplies.

Service: Our worldwide support team is dedicated to providing quality equipment. We walk with you every step of the way to maximize your return on investment:

Made in the USA with a commitment to quality

Warranty: 1 year parts

An ISO 13485:2016 certified company

Standard Features

- Single chamber with windowed door.
- Stainless steel exterior and interior.
- 4 shelves total, 3 adjustable
- Floor-mounted

Capacity

1560-GS
15.6 cubic feet

Configurations (select one)

Door Assembly

- Windowed door (standard)

Size

- 15.6.3 ft³ 15.6 ft³ Built-in
 15.6 ft³ P/T

Door Swing

- Right hinged (standard) Left hinged

Door Lock

- No lock (standard) Manual

Accessories (select all that apply)

Built-in Trim Kit

- Trim kit



Build Your System

Pair your storage cabinet with a warming cabinet. By choosing a left-hinged warming cabinet and a right-hinged storage cabinet, you can create a French-door design and easily transfer supplies from one unit to the next.




Storage cabinet shown with right hinge. EC1350BL shown with optional left hinge

Titan Series Storage Cabinets



Specification



No Clearance Requirements



Model	Exterior (H x W x D)	Interior (H x W x D)	Net Weight
1560-GS	75" (1905mm) x 30" (762mm) x 21.25" (540mm)	60.5" (1537mm) x 24.25" (616mm) x 18" (457mm)	190lbs (86kg)
	Ship Dimensions (L x W x H)*	Ship Weight*	
1560-GS	38" x 35" x 82"	250 lbs	

*Domestic ground shipping information. Contact factory for export weight and dimensions.

CONTACT US

W164 N9221 Water Street | Menomonee Falls, Wisconsin 53051 | U.S.A.
 Phone: 262.251.8356 | 800.862.9276 | Fax: 262.251.7067 | enthermics.com